

## Catalog Number

45016-W

**Tubing Connector Nuts**, low pressure, 1/16 inch, Polypropylene, White with 1/4-28 screw threads, for use with flanged, polymeric tubing. ARE- Applied Research brand. 10/PK.



### Additional Info

These tubing connector nuts are used to connect 1/16th inch (1.6 mm) polymeric tubing in the low pressures zones of HPLC systems and other low to medium pressure, fluidic uses.

- Compatible with 1/4 - 28 UNF fitting ports.
- Useful up to approximately 500 psi or 34.5 bar.
- Maximum Temperature: 100°C.
- Click [HERE](#) for a convenient flanging tool.

Tubing must be flanged, and washers are required to work properly, and sold separately, click [HERE](#) for washer ordering information.

- Maximum Temperature: 100°C.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$42.38**

Printed On 2024/07/20