

## Catalog Number

45099-113

**Polymer Tubing Drill. 1.3mm Polymeric Pin, Replacement for Flanging Tool.** ARE-Applied Research brand. 1 EA.



### Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool.

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$45.68**

Printed On 2024/04/25