

Catalog Number

45099-113-005

Polymer Tubing Drill. 1.3mm Polymeric Pin, Replacement for Flanging Tool. ARE-Applied Research brand. 5/PK.



Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MICROSOLV Flanging Tool.

Price
\$245.20

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$245.20