

Catalog Number

45099-308

Polymer Tubing Drill. 0.8mm Titanium Pin, Replacement for Flanging Tool. ARE-Applied Research brand. Bulk Price.



Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/07/27