

Catalog Number

45099-308-010

Polymer Tubing Drill. 0.8mm Titanium Pin, Replacement for Flanging Tool. ARE-Applied Research brand. 10/PK.



Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool.

Price
\$655.77

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$655.77