

Catalog Number

45099-313

Polymer Tubing Drill. 1.3mm Titanium Pin, Replacement for Flanging Tool. ARE-Applied Research brand. Bulk Pricing.



Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool.

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com