

## Catalog Number

45099-313-005

**Polymer Tubing Drill. 1.3mm Titanium Pin, Replacement for Flanging Tool. ARE-Applied Research brand. 5/PK.**



### Additional Info

This Drill is designed to flange the tips of Polymer (plastic) tubing when used with the MicroSolv Flanging Tool.

**Price**  
**\$344.28**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$344.28**