

## Catalog Number

46000-HP

**Tool Kit for HPLC.** High Pressure. Includes Tools, Tubing, Fittings, Ferrules, Filters and box. ARE-Applied Research brand. 1 EA.



### Additional Info

This convenient kit puts all the parts tools you need to keep an HPLC up and running. A must have for everyone that is responsible for maintenance of the instruments.

Includes:

- 10 Ferrules, PEEK, Double. 1/16"
- 10 Ferrules, SS, Single, 1/16"
- 6 Tubing Connector, One Piece, Finger Tight
- 1 Inline Filter, PEEK complete. 2 um PEEK frit.

**Price**  
**\$899.63**

## Additional Info

- 4 Tubing Connectors, Finger Tight Double Cone Ferrules 10-32
- 4 Tubing Connectors, Hex Head, Long for Double Ferrule 10-32
- 5 Tubing Connectors, SS, Long Head, Rheodyne Type, 10-32
- 5 Tubing Connectors, SS, Extra Long, Rheodyne Type, 10-32
- 2 Solvent Filters, Last Drop, SS Element, 2 um
- 1 Clean Cut Tubing Cutter for Low Pressure Tubing
- 1 File to remove burrs from Tubing.
- 2 Tubing Holders
- 1 Tubing Holder, Knee, 180°
- 1 Tubing Holder, Knee, 90°
- 1 Tubing Lier for cutting metal tubing
- 1 Valve Tool from MicroSolv
- 1 Wrench, 1/4" x 5/16"
- 5 Tubing, PEEK 1/16"OD x 0.18mm ID Striped Yellow
- 5 Tubing, PEEK, 1/16" OD x 0.25mm ID Striped Blue
- 5 Tubing, PTFE, 1/18" x 1/16" ID
- 5 Tubing, SS, 1/16" x 0.25mm ID
- 5 Tubing, SS, 1/16" x 0.50mm ID
- 2 Unions, PEEK, 10-32, 0.3mm Bore, Body Only
- 2 Unions, SS, 10-32, 0.5mm Bore, Body Only

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$899.63**