

## Catalog Number

46000-LP

**Fittings Kit for HPLC.** low pressure. Includes and assortment of flange free fittings and ferrules, in an organizing box. ARE-Applied Research brand. 1 EA.



### Additional Info

This convenient kit puts all the parts you need to keep the low pressure zones of an HPLC up and running. A must have for everyone that is responsible for maintenance of the instruments.

#### Includes:

- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Black 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Blue 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Green 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Red 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", White 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Yellow 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Black 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Blue 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Green 1/4-28

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$617.76**

Printed On 2024/02/29