

Catalog Number

46020-BT

Beta Tool-2. For 1/16 or 1/8" Stainless Steel HPLC Tubing, Blue Grips.
MICROSOLV Brand. 1 EA.



Additional Info

This is a must have tool for every laboratory that works with 1/16" or 1/8" Stainless Steel Tubing. These Tubing Pliers / Tubing Straightener has soft grips and two Center Holes where one is slightly larger than 1/16" and the other is slightly larger than 1/8". This tool is designed to remove kinks or bends in tubing or for making seamless bends when you squeeze the Pliers. This will help to insure a good internal flow path in your HPLC tubing and a void problems due to bent or kinked tubing.

The Beta Tool-2™ features a shorter and rounded nose than most precision pliers and it can reach further into an HPLC Instrument to work with tubing that is already

Price**\$134.71**

Additional Info

installed. For example, working near an injector valve or near some of the tubing that is attached to detectors.

Click [HERE](#) for Product Literature. Click [HERE](#) for a Needle Nose version Beta Tool.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$134.71