


## Catalog Number

46020-FS

**Ferrule Removal & Saver Tool.** For Stainless Steel Ferrules. ARE-  
Applied Research brand. 1 EA.



### Additional Info

When a Stainless Steel Ferrule swages onto tubing, it is very difficult to remove it without damaging the Ferrule or the Tubing.

This tool is designed to remove CPI type Stainless Steel Ferrules from swaged HPLC Tubing, easily and without damaging the tubing.

*Helps avoid cutting Tubing and making tubing too short or cutting off a perfect cut tip.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$175.50**

Printed On 2024/05/04