

## Catalog Number

46040-02

**Clean Cut90 Tubing Cutter**, for PEEK Tubing, with Storage Case & 1 Replacement Blade. ARE-Applied Research brand. 1 EA.



### Additional Info

This is one of our most popular tubing cutters since it is easy to operate properly and produces a very good, nearly 90° cut. When cut properly, the PEEK tubing can be used in critical HPLC applications. While mainly used for PEEK tubing, other polymeric tubing such as PTFE can be cut as well.

The Clean-Cut90™ can cut 5 different tubing outer diameters such as 1/4 inch, 0.15 inch, 1/8 inch, 0.105 inch and 1/16 inch. The cutting blade has a "flat side" and a "back side". An indicating arrow shows the side of the tubing that will have the near 90, perpendicular cut for use on your instrument which is the cutting side. Click [HERE](#) for use instructions.

*A plastic storage case and 1 replacement blade is included.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$28.89**

Printed On 2024/05/04