


## Catalog Number

46100-01T-005

**Tweezers**, Polypropylene, Blue. ARE-Applied Research brand. 5/PK.



### Additional Info

These lightweight, plastic Tweezers are handy for picking light weight objects that should not be touched by human hands to avoid contamination.

Useful for Membrane Filters, Sparging Stones, Weighing Paper and many other items around a Laboratory.

**Price**  
**\$12.78**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$12.78**