



Catalog Number

46200-01F

HPLC File, for Stainless Steel Tubing. ARE-Applied Research brand.
Bulk Price.



Additional Info

This file is used when cutting or breaking Stainless Steel Tubing for use in HPLC. The proper finish on the file helps to remove burrs and to smooth tubing to achieve a satisfactory tip for a proper connection.



Additional Info

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com