

Catalog Number

46200-01F

HPLC File, for Stainless Steel Tubing. ARE-Applied Research brand. 1 EA. M.O.Q. 11 Pieces.



Additional Info

This file is used when cutting or breaking Stainless Steel Tubing for use in HPLC. The proper finish on the file helps to remove burrs and to smooth tubing to achieve a satisfactory tip for a proper connection.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$17.34

Printed On 2024/02/25