

## Catalog Number

46200-01F-005

**HPLC File**, for Stainless Steel Tubing. ARE-Applied Research brand.  
5/PK.



### Additional Info

This file is used when cutting or breaking Stainless Steel Tubing for use in HPLC. The proper finish on the file helps to remove burrs and to smooth tubing to achieve a satisfactory tip for a proper connection.

Easy order pack of 5. Click [HERE](#) for pack of 10. Click [HERE](#) if you need a larger quantity.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$91.00**

Printed On 2024/05/18