

## Catalog Number

46552-180

**Tubing Holder, Knee, 180 Degree Turn.** Blue. ARE-Applied Research brand. 10/PK.



### Additional Info

This Tubing Holder, will keep your tubing "in place" when making a 180° turn, which is not always easy. Avoid clutter and keep your system optimized.

These tubing knees will give you the right bend of your PEEK tubing, preventing close diameters, which do have a negative effect on the mechanical stability of your PEEK tubing and can cause poor internal dimensions.

Best performance when used with 1/16th" PEEK Tubing but other tubing can be used as well.

**Price**  
**\$29.54**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$29.54**