

Catalog Number

46552-90

Tubing Holder, Knee, 90 Degree Turn. Blue. ARE-Applied Research brand. 10/PK.



Additional Info

This Tubing Holder, will keep your tubing "in place" when making a 90° turn, which is not easy. Avoid clutter and keep your system optimized.

These tubing knees will give you the right bend of your PEEK tubing, preventing close diameters, which do have a negative effect on the mechanical stability of your PEEK tubing and can cause poor internal dimensions.

Best performance when used with 1/16th" PEEK Tubing but other tubing can be used as well.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$28.89

Printed On 2024/05/04