


## Catalog Number

46600-01

**Tubing Holder, Clamp**, for 1/8" or 1/16" OD Tubing. Blue. ARE-Applied Research brand. 10/PK.



### Additional Info

This Tubing Holder will hold 1/16" or 1/8" OD Tubing in place on a beaker, bottle or any thin walled vessel. The spring loaded holder is excellent while draining or siphoning or filling with tubing.

Equipped with 2 holes, one hole accommodates 1/16" and the other accommodates 1/8" tubing.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$48.15**

Printed On 2024/02/25