


# Catalog Number

47305-05

**Tubing, PEEK, 0.005 inch (0.13mm) ID, 1/16th inch (1.6 mm) OD, HPLC grade, red stripe, 5 meter roll, ARE-Applied Research brand. 1 EA.**



## Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The red stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.735 mm
- Color code: red stripe
- Pressure rating: 6100 psi (420.6 bar)

**Price**  
**\$96.30**

## Additional Info

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$96.30**