


Catalog Number

47305-05-10M

Tubing, PEEK, 0.005 inch (0.13mm) ID, 1/16th inch (1.6 mm) OD,
HPLC grade, red stripe, 10 meter roll, ARE-Applied Research brand. 1
EA.



Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The red stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.735 mm
- Color code: red stripe
- Pressure rating: 6100 psi (420.6 bar)

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$192.60

Printed On 2024/05/07