


## Catalog Number

47305-05-1C

**Tubing, PEEK, 0.005 inch (0.13 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, red stripe, Your custom continuous length tubing, priced per  
meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this HPLC grade tubing that is extruded from virgin PEEK material without any manufacturing additives. The red stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.735 mm
- Color code: red stripe
- Pressure rating: 6100 psi (420.6 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/18