




Catalog Number

47305-10-3M

Tubing, PEEK, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD, HPLC grade, blue stripe, 3 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The blue stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.675 mm
- Color code: blue stripe
- Pressure rating: 5600 psi (386.1 bar)

Price
\$82.06

Additional Info

Other lengths available:

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

Custom cuts also available. Please contact [Customer Service](#).

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$82.06