


Catalog Number

47305-13-10

**Tubing, PEEK, Pre-Cut, 0.005 inch (0.13 mm) ID, 1/16th inch (1.6 mm) ID, Ultra-T grade, red intermittent stripes, 10 cm long (+/- 1 mm). ARE-Applied Research brand. 5 PK.**



### Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut PEEK Tubes are the perfect solution.

This Ultra-T™ grade (*Premium*), PEEK tubing is extruded from virgin PEEK material without adding any additives. The broken red stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color. The internal and outer dimensional tolerances of this tubing is +/- 0.025mm.

- Color Code: Red, Stripes for 0.005" (0.13mm) ID.

**Price**  
**\$38.53**

## Additional Info

- Wall Thickness: 0.735 mm
- Pressure rating: 6100 psi.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$38.53**