


Catalog Number

47305-25-10

Tubing, PEEK, Pre-Cut, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD, Ultra-T grade, blue intermittent stripes, 10 cm (+/- 1 mm). ARE- Applied Research brand. 5 PK.



Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut PEEK Tubes are the perfect solution.

This Ultra-T™ grade (*Premium*), PEEK tubing is extruded from virgin PEEK material without adding any additives. The broken yellow stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color. The internal and outer dimensional tolerances of this tubing is +/- 0.025mm.

- Color Code: Blue, Stripes for 0.010" (0.25mm) ID.
- Wall Thickness: 0.675 mm

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$38.53

Printed On 2024/08/06

Additional Info

- Pressure rating: 5600 psi.

NOTE: *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$38.53

Printed On 2024/08/06