


Catalog Number

47305-30

Tubing, PEEK, 0.030 inch (0.75 mm) ID, 1/16th inch (1.6 mm) OD, HPLC grade, green stripe, 5 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The green stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.425 mm
- Color code: green stripe
- Pressure rating: 3500 psi (241.3 bar)

Other lengths available:

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$96.30

Printed On 2024/08/06

Additional Info

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

NOTE: This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$96.30

Printed On 2024/08/06