


Catalog Number

47305-40-10M

Tubing, PEEK, 0.040 inch (1.0 mm) ID, 1/16th inch (1.6 mm) OD, HPLC grade, grey stripe, 10 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The grey stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.30 mm
- Color code: grey stripe
- Pressure rating: 2400 psi (165.5 bar)

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$234.00

Printed On 2024/05/18