

## Catalog Number

47305-55-1M

**Tubing, PEEK, 0.055 inch (1.4 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, black stripe, 1 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The black stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.10 mm
- Color code: black stripe
- Pressure rating: 750 psi (51.7 bar)
- *Caution Note: This tubing has very thin walls and is easily damaged. Handle with care to prevent kinking.*

Other lengths available:

- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$27.15**

Printed On 2024/04/25