

Catalog Number

47305-ST-02-1M

Tubing, PEEK, 0.0025 inch (0.064 mm) ID, 1/16th inch (1.6 mm) OD, Ultra-T grade, natural color, 1 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. A natural PEEK material has no color leachables.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.768 mm
- Color code: Natural (no stripe)
- Pressure rating: 6300 psi (434.4 bar)

Other lengths available:

Price
\$34.83

Additional Info

- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

NOTE: This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$34.83