


## Catalog Number

47305-ST-02-C

**Tubing, PEEK, 0.0025 inch (0.064 mm) ID, 1/16th inch (1.6 mm) OD, Ultra-T grade, natural color, 25 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. A natural PEEK material has no color leachables.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.768 mm
- Color code: Natural (no stripe)
- Pressure rating: 6300 psi (434.4 bar)
- Supplied on a convenient spool.

**Price**  
**\$787.73**

## Additional Info

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [Custom lengths](#)

**NOTE:** *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$787.73**