


## Catalog Number

47305-ST-05-1C

**Tubing, PEEK, 0.005 inch (0.13 mm) ID, 1/16th inch (1.6 mm) OD,**  
Ultra-T grade, natural color, Your custom continuous length tubing, priced  
per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this Ultra-T™ grade (*premium*) tubing that is extruded from virgin PEEK material without any manufacturing additives.

The broken red stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.735 mm
- Color code: Red stripes
- Pressure rating: 6100 psi (420.6 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/04/28