


## Catalog Number

47305-ST-07-1C

**Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/16th inch (1.6 mm) OD,** Ultra-T grade, yellow intermittent stripes, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this Ultra-T™ grade (*premium*) tubing that is extruded from virgin PEEK material without any manufacturing additives. The broken yellow stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.71 mm
- Color code: yellow stripes
- Pressure rating: 5800 psi (399.9 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/07