


Catalog Number

47305-ST-10-1C

Tubing, PEEK, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD,
Ultra-T grade, blue intermittent stripes, Your custom continuous length
tubing, priced per meter. ARE-Applied Research brand.



Additional Info

You determine the length of this Ultra-T™ grade (*premium*) tubing that is extruded from virgin PEEK material without any manufacturing additives. The broken blue stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.675 mm
- Color code: blue stripes
- Pressure rating: 5600 psi (386.1 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/03