

Catalog Number

47310-05

Tubing, Capillary, PEEK, 75 um ID, 360 um OD, natural color, 10 Meter Roll. ARE-Applied Research brand. 1 EA.



Additional Info

This small OD, PEEK, capillary tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.157 um
- Color code: natural
- Pressure rating: 5200 psi (358.5 bar)

Other lengths available:

Price
\$295.40

Additional Info

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)
- [25 meter roll](#)
- [Custom lengths](#)

NOTE: This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$295.40