


## Catalog Number

47310-100-1C

**Tubing, PEEK, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD,**  
HPLC grade, natural color, Your custom continuous length tubing, priced  
per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this wide OD, PEEK, tubing. It is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.80 mm
- Color code: natural
- Pressure rating: 3200 psi (220.6 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/08/06

## Additional Info

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

**NOTE:** *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/08/06