




Catalog Number

47310-100-5M

Tubing, PEEK, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD, HPLC grade, natural color, 5 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.80 mm
- Color code: natural
- Pressure rating: 3200 psi (220.6 bar)

Other lengths available:

Price

\$194.95



Additional Info

- [1 meter roll](#)
- [3 meter roll](#)
- [10 meter roll](#)
- [25 meter roll](#)

Custom cuts also available. Please contact [Customer Service](#).

NOTE: *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

MICROSOLV

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$194.95