

Catalog Number

47310-13-1C

Tubing, PEEK, 0.005 inch (0.13 mm) ID x 1/32nd inch (0.79 mm) OD, general grade, natural, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.



Additional Info

You determine the length of this general grade tubing that is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.33 mm
- Color code: Natural
- Pressure rating: 5500 psi (380 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/07/27