

## Catalog Number

47310-13-1C

**Tubing, PEEK, 0.005 inch (0.13 mm) ID x 1/32nd inch (0.79 mm) OD,** general grade, natural, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this general grade tubing that is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.33 mm
- Color code: Natural
- Pressure rating: 5500 psi (380 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/24