

## Catalog Number

## 47310-18-C

**Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/32nd inch (0.79 mm) OD, general grade, natural, 25 meter roll, ARE-Applied Research brand. 1 EA.**



### Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.305 mm
- Color code: Natural
- Pressure rating: 5100 psi (350 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$866.70**

Printed On 2024/05/18