

Catalog Number

47310-200-C

Tubing, PEEK, 0.080 inch (2.0 mm) ID, 1/8th inch (3.2 mm) OD, HPLC grade, natural color, 25 meter roll, ARE-Applied Research brand. 1 EA.



Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.60 mm
- Color code: natural
- Pressure rating: 2400 psi (165.5 bar)
- Supplied on a convenient spool.

Price
\$625.95

Additional Info

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)
- [Custom lengths](#)

NOTE: *This tubing is cut with great care to be even but is cut for length but not as a precision cut required by some HPLC situations. Please inquire for more information if you would like your tubing cut to a length ready for use in low dead volume HPLC fittings.*

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$625.95