


Catalog Number

47310-300-C

Tubing, PEEK, 1/8th inch (3.2 mm) ID, 1/4th inch (6.35 mm) OD,
HPLC grade, natural color, 25 meter roll, ARE-Applied Research brand. 1
EA.



Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 1.575 mm
- Color code: natural
- Pressure rating: 3300 psi (227.5 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$4,526.10

Printed On 2024/05/18