


## Catalog Number

47320-30

**Tubing connector fittings**, high pressure, two pieces, 1/16 inch, PEEK, natural, large head, standard length and 10-32 screw threads. Use with double cone ferrules. Finger tighten, ARE-Applied Research brand. Bulk



**Bulk Pricing**

**MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)