


Catalog Number

47320-31-005

Tubing connector fittings, high pressure, two pieces, 1/16 inch, PEEK, natural, large head, standard length and 10-32 screw threads. Use with double cone ferrules. Finger tighten, ARE-Applied Research brand. 5/PK.



Pack of 5

Price
\$45.68

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price**\$45.68**