

Catalog Number

47452-04

Union, for use with low pressure. Barbed connection for 1/16 inch tubing. Polypropylene. ARE-Applied Research brand. Bulk Price.



Additional Info

This union connects 1/6th inch (1.6 mm) ID tubing when slid over the barbed fittings. *Bulk packaging.*

- Maximum Temperature: 100°C
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for economical packs of 10.

For a tubing adapter to connect 1/16th inch (1.6 mm) ID tubing to 1/8th inch (3.2 mm) ID tubing, Click [HERE](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/17