

Catalog Number

47452-11-010

Adapter, for use with low pressure HPLC. Female luer to 1/4-28 male connection. Polypropylene. ARE-Applied Research brand. 10/PK.



Additional Info

This adapter has a female Luer Slip adapter on one side and a 1/4-28 threaded male connector on the other side. Bulk packaging.

- Maximum Temperature: 100C
- Maximum Pressure: 1,000 psig.

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

Price
\$43.21

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$43.21