

## Catalog Number

47452-12-010

**Adapter**, for use with low pressure. Female luer lock to barbed connection for 1/8 inch tubing. Polypropylene. ARE-Applied Research brand. 10/PK.



### Additional Info

This adapter has a female Luer Lock port on one side and a hose barb connector on the other side for 1/8th inch (3.2 mm) ID tubing.

This adapter is excellent when connecting 1/8th inch tubing to a syringe.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

**Price**  
**\$42.26**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$42.26**