

Catalog Number

47452-13-010

Adapter, for use with low pressure HPLC. Female luer to 10/32 male connection. Polypropylene. ARE-Applied Research brand. 10/PK.



Additional Info

This adapter has a female Luer Lock port on one side and a 10-32 threaded male connector on the other side.

This adapter is excellent when using a syringe and you want to connect to an instrument port with a 10-32 thread.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

Price
\$43.21

Additional Info

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$43.21