

Catalog Number

47452-32-010

Adapter, for use with low pressure. Male luer to barbed connection for 1/8 inch tubing. Polypropylene. ARE-Applied Research brand. 10/PK.



Additional Info

This adapter has a male Luer Lock connection on one side and a hose barb male connection on the other side that accommodates 1/8th inch (3.2 mm) ID tubing.

This adapter is excellent when using a syringe filter or other device that has a female Luer Lock or Luer Slip port and you want to connect with 1/8th inch tubing

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

Price
\$42.26

Additional Info

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$42.26