



# Catalog Number

47453-01-010

Adapter, for use with low pressure HPLC. Female luer thread style with 5/16 inch Hex to 1/4-28 UNF bottom sealing thread. White nylon ARE- Applied Research brand. 10/PK.



## Additional Info

This adapter has Luer fitting and 1/4-28 thread for sealing. ID tubing on both sides.  
*Bulk packaging.*

*This adapter is great when tubing needs different connectors.*

- Maximum Temperature: 100°C
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for Bulk Pricing.



## Additional Info

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

### **MICROSOLV**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)