

## Catalog Number

47453-02-010

Coupler, for use with low pressure HPLC. Male Luer Coupler. White nylon ARE-Applied Research brand. 10/PK.

**Additional Info**

A high-quality, white nylon connector designed for secure and reliable connections in medical and laboratory applications. This coupler is compatible with separate rotating lock rings, providing flexibility and ease of use in various setups.

**Key Features:**

- **Material:** Made from durable white nylon, ensuring longevity and resistance to wear and tear.
- **Compatibility:** Can be used with separate rotating lock rings for enhanced

**Price**  
**\$13.19**

## Additional Info

connection stability.

- **Design:** Features a standardized 6% luer taper, ensuring compatibility with other luer fittings.
- **Applications:** Ideal for use in laboratory instruments, and fluid management systems where secure, leak-proof connections are critical.
- **Tubing Connections:** Suitable for tubing with inner diameters ranging from 1/16 inch (1.6 mm) to 1/4 inch (6.4 mm), providing versatility in various applications.

## Benefits:

- **Leak-Proof:** The precise design ensures a secure friction seal, preventing leaks in both gas and fluid connections.
- **Versatile:** Compatible with a wide range of tubing materials and sizes, making it suitable for laboratory applications.
- **Reliable:** Manufactured to high standards, ensuring consistent performance and reliability in critical applications.

## Usage:

- **Laboratory Instruments:** Ideal for use in analytical instruments, fluid transfer systems, and other laboratory equipment.

## MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

## Price

**\$13.19**